

Contec® NeutraKlean

Irradiated neutral detergent with purified water

Contec® NeutraKlean is a gentle and non-foaming detergent designed for the cleaning of life science cleanrooms.

With a pH between 6.0 and 8.0 after irradiation, Contec NeutraKlean is ideal for carrying out the low-level cleaning required after maintenance, and to remove disinfectant residues and general soil. The use of a surfactant/detergent helps to break the binding of particles and dirt to a surface and allows them to be more easily removed.

The product is 0.2 µm filtered and filled under ISO Class 5 airflow and packaged in an ISO Class cleanroom. Ready to use, the product is available in a 1 L trigger spray with a “bag-in-bottle” system to protect the contents during use. A non-sterile version is available in a 5 L capped container.



Both the sterile and filtered product are provided double bagged to aid the transfer of the product into a pharmaceutical cleanroom. The linear tear bags are easy to open even when wearing gloves.



Features	Benefits
Neutral pH	<ul style="list-style-type: none"> • Suitable for use in all areas of cleanroom and gentle on surfaces
Filtered to 0.2 microns under ISO Class 5 airflow in a ISO Class 7 cleanroom	<ul style="list-style-type: none"> • Ensures the product is free from contamination and particulates
Low foaming detergent action	<ul style="list-style-type: none"> • Ideal for use on floors, walls, and benches
Low level of surfactant	<ul style="list-style-type: none"> • Low residue which can be easily removed
Trigger spray can be set to stream or mist	<ul style="list-style-type: none"> • Large droplet size reduces the risk of inhalation and provides good surface coverage
Trigger spray and “bag-in-bottle” protected system	<ul style="list-style-type: none"> • Prevents bottle contents from becoming contaminated during use
Double bagged packed in linear tear packaging	<ul style="list-style-type: none"> • Each bag is easy to open even when wearing gloves • Facilitates transfer of cleaning products into isolators and RABS

Full documentation available upon request.


Part No.	Description		Size	Packaging
SBT103NK	Contec Irradiated NeutraKlean Trigger spray		34 oz. (1 L)	6/case
FBC503NK	Contec NeutraKlean		1.3 gallon. (5 L)	2/case


Product Information					
pH after irradiation	• 6.0 - 8.0				
Packaging materials	• Outer bags (OB1, OB2), low density polyethylene (LDPE)  Case (CS), corrugated fiberboard (PAP) 				
Environment	• ISO 5-8 Grade A/B for sterile, C/D for nonsterile				
Shelf life	<table border="0"> <tr> <td>Unopened</td> <td>• 2 years</td> </tr> <tr> <td>Opened</td> <td>• 6 months</td> </tr> </table>	Unopened	• 2 years	Opened	• 6 months
Unopened	• 2 years				
Opened	• 6 months				


Recycle Symbols

- PET 
- HDPE 
- LDPE 
- PP 
- PAP 

Contec® Disinfectant and Solution Rotation

DAILY	Use everyday on all surfaces.	Disinfectants CyQuanol™ Sterile CyQuanol™ PREempt® Plus PREempt® Concentrate	
--------------	-------------------------------	---	--

ROUTINELY	Use based on room environment monitoring.	Sporicide PeridoxRTU® Sterile PeridoxRTU Peridox® Concentrate	
------------------	---	---	---

PERIODICALLY	Use to remove built up residue or when needed.	Detergent NeutraKlean Sterile NeutraKlean	
---------------------	--	--	---

RINSE		
Use after cleaning and disinfecting.	Rinse IPA HydroKlean WFI Water	