

# Contec® NeutraKlean

## Irradiated neutral detergent with purified water

Contec® NeutraKlean is a gentle and non-foaming detergent designed for the cleaning of life science cleanrooms.

With a pH between 6.5 and 8.0 after irradiation, Contec NeutraKlean is ideal for carrying out the low-level cleaning required after maintenance, and to remove disinfectant residues and general soil. The use of a surfactant/detergent helps to break the binding of particles and dirt to a surface and allows them to be more easily removed.

The product is 0.2 µm filtered and filled under ISO Class 5 airflow and packaged in an ISO Class cleanroom. Ready to use, the product is available in a 1 L trigger spray with a “bag-in-bottle” system to protect the contents during use. A non-sterile version is available in a 5 L capped container.



Both the sterile and filtered product are provided double bagged to aid the transfer of the product into a pharmaceutical cleanroom. The linear tear bags are easy to open even when wearing gloves.



Features	Benefits
Neutral pH	<ul style="list-style-type: none"> <li>• Suitable for use in all areas of cleanroom and gentle on surfaces</li> </ul>
Filtered to 0.2 microns under ISO Class 5 airflow in a ISO Class 7 cleanroom	<ul style="list-style-type: none"> <li>• Ensures the product is free from contamination and particulates</li> </ul>
Low foaming detergent action	<ul style="list-style-type: none"> <li>• Ideal for use on floors, walls, and benches</li> </ul>
Low level of surfactant	<ul style="list-style-type: none"> <li>• Low residue which can be easily removed</li> </ul>
Trigger spray can be set to stream or mist	<ul style="list-style-type: none"> <li>• Large droplet size reduces the risk of inhalation and provides good surface coverage</li> </ul>
Trigger spray and “bag-in-bottle” protected system	<ul style="list-style-type: none"> <li>• Prevents bottle contents from becoming contaminated during use</li> </ul>
Double bagged packed in linear tear packaging	<ul style="list-style-type: none"> <li>• Each bag is easy to open even when wearing gloves</li> <li>• Facilitates transfer of cleaning products into isolators and RABS</li> </ul>

Full documentation available upon request.


Part No.	Description		Size	Packaging
SBT103NK	Contec Irradiated NeutraKlean Trigger spray		34 oz. (1 L)	6/case
FBC503NK	Contec NeutraKlean		1.3 gallon. (5 L)	2/case


Product Information					
pH after irradiation	• 6.5 - 8.0				
Packaging materials	• Outer bags (OB1, OB2), low density polyethylene (LDPE)  Case (CS), corrugated fiberboard (PAP) 				
Environment	• ISO 5-8 Grade A/B for sterile, C/D for nonsterile				
Shelf life	<table border="0"> <tr> <td>Unopened</td> <td>• 2 years</td> </tr> <tr> <td>Opened</td> <td>• 6 months</td> </tr> </table>	Unopened	• 2 years	Opened	• 6 months
Unopened	• 2 years				
Opened	• 6 months				


Recycle Symbols

- PET 
- HDPE 
- LDPE 
- PP 
- PAP 

**Contec® Disinfectant and Solution Rotation**

<b>DAILY</b>	Use everyday on all surfaces.	<b>Disinfectants</b> <ul style="list-style-type: none"> <li>• CyQuano™</li> <li>• Sterile CyQuano™</li> <li>• PREempt™ RTU</li> <li>• PREempt™ Concentrate</li> </ul> 
--------------	-------------------------------	--

<b>ROUTINELY</b>	Use based on room environment monitoring.	<b>Sporicide</b> <ul style="list-style-type: none"> <li>• PeridoxRTU®</li> <li>• Sterile PeridoxRTU®</li> <li>• Peridox® Concentrate</li> </ul> 
------------------	---	---

<b>PERIODICALLY</b>	Use to remove built up residue or when needed.	<b>Detergent</b> <ul style="list-style-type: none"> <li>• NeutraKlean</li> <li>• Sterile NeutraKlean</li> </ul> 
---------------------	--	---

<b>RINSE</b>	
Use after cleaning and disinfecting.	<b>Rinse</b> <ul style="list-style-type: none"> <li>• IPA</li> <li>• HydroKlean</li> <li>• WFI</li> <li>• Water</li> </ul> 