

# Contec® Irradiated 70% isopropyl alcohol LE



## Low endotoxin irradiated isopropanol with water for injection



Contec® Low Endotoxin Irradiated Alcohol contains 70% by volume USP Grade isopropyl alcohol (IPA) and 30% Water for Injection (WFI). The alcohol blend is 0.2 micron filtered, filled, and bagged in an ISO Class 7 cleanroom. This clean manufacture coupled with water for injection means the alcohol blend is guaranteed and certified to have an endotoxin level < 0.25 EU/ml, making it suitable for product contact surfaces. The bottles are provided with the triggers installed and gamma irradiated at no less than 25kGy. Each lot is sterility tested and certified.

Available as a 17 oz and 34 oz trigger spray bottle, the trigger sprayers are fitted with a protected system, which ensure sterility throughout use. Each bottle has a lot code and expiration date for easy record keeping. The product has a 3 year unopened shelf-life from date of manufacture and a 6 month in-use shelf-life.



Features	Benefits
70% isopropyl alcohol and 30% Water for Injection	<ul style="list-style-type: none"> <li>• No residue for product contact areas</li> </ul>
Guaranteed endotoxin level of less than 0.25 EU/ml	<ul style="list-style-type: none"> <li>• Ideal for use in product contact areas</li> </ul>
Double bagged in linear tear packaging	<ul style="list-style-type: none"> <li>• Easy to open even when wearing gloves</li> <li>• Facilitates transfer into isolators and RABS</li> </ul>
Gamma irradiated at no less than 25kGy and lot tested for sterility	<ul style="list-style-type: none"> <li>• Suitable for use in the most critical environments</li> </ul>
Trigger spray in “bag-in-bottle” protected system	<ul style="list-style-type: none"> <li>• Eliminates the risk of the alcohol being contaminated when installing the trigger</li> <li>• Prevents bottle contents from being contaminated during use</li> </ul>
Lot numbered bottles	<ul style="list-style-type: none"> <li>• Assures complete product traceability</li> </ul>
Trigger sprayer can be set to stream or mist	<ul style="list-style-type: none"> <li>• Large droplet size reduces the risk of inhalation and provides good surface coverage</li> </ul>

Part No.	Description		Size	Packaging
SBT167030LE	Contec Irradiated 70% IPA LE Sterile 70% isopropyl alcohol in water for injection trigger spray		17 oz. (0.5 L)	8/case
SBT347030LE	Contec Irradiated 70% IPA LE Sterile 70% isopropyl alcohol in water for injection trigger spray		34 oz. (1 L)	6/case

Product Information	
Appearance	<ul style="list-style-type: none"> <li>• Clear, colorless liquid</li> </ul>
Odor	<ul style="list-style-type: none"> <li>• Characteristic alcohol</li> </ul>
% IPA Contents	<p>70%</p> <ul style="list-style-type: none"> <li>• 68.0% - 72.0%</li> </ul>
Packaging materials	<ul style="list-style-type: none"> <li>• Outer bags (OB1, OB2), low density polyethylene (LDPE) </li> <li>Case (CS), corrugated fiberboard (PAP) </li> </ul>
Environment	<ul style="list-style-type: none"> <li>• ISO 3-8 Grade A/B</li> </ul>
Shelf life	<ul style="list-style-type: none"> <li>Unopened • 3 years</li> <li>Opened • 6 months</li> </ul>

Recycle Symbols



Packaging	EA/OB1	OB1/OB2	OB2/CS	EA/CS
SBT167030LE	1	1	8	8
SBT347030LE	1	1	6	6

EA = each; OB = outer bag; CS = case