

PREempt® Plus

One-step disinfectant cleaner utilizing Accelerated Hydrogen Peroxide®

PREempt® Plus is a formulation of Accelerated Hydrogen Peroxide® (AHP®) that features faster kill times and lower residue. PREempt Plus is a Ready-to-Use, one-step broad-spectrum cleaner and disinfectant designed for cleaning and disinfecting hard, non-porous surfaces.

PREempt Plus kills viruses, bacteria, fungi, and TB in 1 minute or less.*

PREempt Plus has lower residue than the original PREempt RTU making it more suitable for critical environments. It contains no volatile organic compounds (VOCs) and carries no hazard classifications.

PREempt Plus utilizes Accelerated Hydrogen Peroxide (AHP), a patented synergistic blend of commonly used ingredients, that when combined with low levels of hydrogen peroxide dramatically increase its efficacy and cleaning performance.

PREempt Plus products are ideal for cleaning and disinfecting environmental surfaces in cleanrooms and laboratory areas including workstations, fume hoods, equipment and other hard non-porous environmental surfaces.



Features	Benefits
Faster kill times	<ul style="list-style-type: none"> • Easy to maintain compliance and appropriate dwell time
Less residue	<ul style="list-style-type: none"> • Less than half the residue of traditional PREempt RTU • More suitable for critical environments
Contains low levels of hydrogen peroxide	<ul style="list-style-type: none"> • Breaks down into water and oxygen so low residue • Nonhazardous so allows cost effective disposal
Lowest EPA toxicity rating	<ul style="list-style-type: none"> • Odor-free for better operator acceptability • Contains no VOCs • Requires no PPE
Compatible with a wide range of materials	<ul style="list-style-type: none"> • Reduces corrosion and wear on surfaces and equipment

Manufactured by:



Engineering Revolutionary Disinfectants for the War Against Microbes



Accelerated Hydrogen Peroxide®

PDSD017 | 120622

Copyright © 2022 Contec, Inc. All rights reserved.

Contact us:
 info@contecinc.com
 contecinc.com

Contec, Inc.
 Spartanburg, SC 29303 USA
 tel: +1-864-503-8333



Part No.	Description	Size	Packaging
2B101	PREempt Plus Includes 12 flip-top caps	32 oz. (0.9 L)	12/case
2B105	PREempt Plus Capped	1 gallon (3.8 L)	4/case
2B221	PREempt Plus Wipes Meltblown polypropylene wipes	6" x 7" (152 x 178 mm)	12/case
ACCTRIGR	Non-aerosol trigger sprayer for use with 32 oz. bottle	-	12/case

Product Application:

PREempt Plus is a one-step cleaner/disinfectant used to treat hard, non-porous surfaces. If area is visibly soiled, precleaning may be required. Apply PREempt Plus to surfaces using a cleanroom wipe or mop. Surfaces must remain wet for the duration of the dwell time. PREempt Plus can be rinsed with Deionized Water (DI) or Water for Injection (WFI), 70% Isopropyl Alcohol (IPA)/Purified Water (PW) or 70% IPA/WFI blends, or dry wipe after application and achieving the required contact time. Always follow site Standard Operating Procedures (SOPs).

It is a violation to use this product in a manner inconsistent with its labelling.

Specification

- Contains 0.5% Hydrogen Peroxide
- 3-year shelf life (opened and unopened)
- Relevant certificates of analysis available for every batch
- Additional supporting documentation is available on request

Technical Data			
Test	Organism	Species	Time
EPA Sanitizer Test for Non-Food Contact Surfaces	Bacteria	<i>E.coli</i> <i>P. aeruginosa</i> <i>S. aureus</i>	10 secs
AOAC Use Dilution Method	Bacteria	<i>E.coli</i> <i>P. aeruginosa</i> <i>S. aureus</i>	1 min
EPA Quantitative Tuberculocidal Activity Test Method	TB	<i>Mycobacterium bovis</i> (BCG)	1 min
Standard Test Method for Efficacy of Virucidal Agents	Viruses	<i>Poliovirus Type 1</i> <i>Adenovirus Type 8</i> <i>Rhinovirus Type 37</i>	30 secs*
AOAC Fungicidal Activity Disinfectants Model	Fungi	<i>Candida albicans</i>	1 min

*Animal viruses may require contact times of 30 seconds - 3 minutes. Please refer to the product label for more information.