

CONSTIX® Sealed Polyester Swab SP-6M


Double layer sealed knitted microfiber swab



CONSTIX® SP-6M sealed polyester swab consists of 100% knitted polyester microfiber fabric that is thermally bonded to a long flex-tip polypropylene handle. Constructed without adhesives or binders, the swab head is made of cleanroom laundered, double layer knitted microfiber polyester materia and is compatible with most solvents, such as alcohols and ketones. The flexible tip of the SP-6M swab provides a softer application of the swab head to surfaces. The swab is durable, absorbent and features extremely low particles and residue levels.



Features	Benefits
Sealed polyester fabric	<ul style="list-style-type: none"> • Provides enhanced sorbency, particle entrapment, and removal • Compatible with most common solvents, such as alcohols and ketones
Flexible tip	<ul style="list-style-type: none"> • Provides a softer application on surfaces

Part No.	Description	Size	Packaging
SP-6M	CONSTIX Sealed polyester swab	Head dimension: 0.27" x 0.64" (6.9 x 16.3 mm) Handle diameter: 0.12" (3.0 mm) Total length: 6.4" (162.6 mm)	500/case



Product Information	
Head material	<ul style="list-style-type: none"> Knitted microfiber polyester
Handle material	<ul style="list-style-type: none"> Polypropylene
Construction	<ul style="list-style-type: none"> Thermally sealed and attached
Packaging materials	<ul style="list-style-type: none"> Outer bags (BG): low density polyethylene (LDPE)  Case (CS): corrugated fiberboard (PAP) 

Recycle Symbols



Packaging	EA/BG	BG/CS	EA/CS
SP-6M	100	5	500

EA = each; BG = bag; CS = case