



# PROSAT Sterile Low Endotoxin Wipes

## Polynit Heatseal wipes with 70% IPA & Water for Injection

PROSAT Sterile Low Endotoxin Wipes are 100% knitted polyester Polynit Heatseal wipes presaturated with a blend of 70% IPA and Water for Injection (EP). Each batch of wipes has a guaranteed endotoxin level of less than <1EU per wipe.

The sealed edge wipe produces a very low level of particles and fibres and is ideal for use on product contact surfaces. The hanging slot in the pouch makes it easy to use in isolators and RABS.

Each pouch only contains a small number of wipes eliminating any waste at the end of a session. Double or triple bagged and validated sterile, the pouches are suitable for use in the highest grades of pharmaceutical cleanroom.



FEATURES	BENEFITS
<b>Cleanroom laundered 100% knitted polyester wipes</b>	<ul style="list-style-type: none"> <li>• Very low levels of particles and fibres</li> <li>• Good abrasion and chemical resistance</li> </ul>
<b>Laser sealed edges</b>	<ul style="list-style-type: none"> <li>• Reduces fiber generation from the edge of the wipe</li> </ul>
<b>Presaturated with 70% IPA/WFI (EP)</b>	<ul style="list-style-type: none"> <li>• Controlled saturation levels greatly reduce solvent usage, VOC's &amp; waste</li> </ul>
<b>Guaranteed endotoxin levels of &lt;1EU/wipe</b>	<ul style="list-style-type: none"> <li>• Suitable for use on product contact surfaces</li> </ul>
<b>Validated sterile via gamma irradiation to a 10<sup>-6</sup> Sterility Assurance Level</b>	<ul style="list-style-type: none"> <li>• Suitable for use in Grade A and B cleanrooms</li> </ul>
<b>Small number of wipes per resealable pouch</b>	<ul style="list-style-type: none"> <li>• Eliminates wastage.</li> <li>• Pouch can be easily resealed throughout use, maintaining saturation level</li> </ul>
<b>Double or triple bagged in linear tear packaging</b>	<ul style="list-style-type: none"> <li>• Each bag is easy to open even when wearing gloves</li> <li>• Facilitates cleaning products entering isolators and RABS</li> </ul>

Part No.	Description	Size	Quantity
PSWE0001	PROSAT Sterile Polynit Heatseal LE Wipes Presaturated with 70% IPA and 30% WFI (EP), triple bagged	230 x 230mm	10 wipes per pouch; 55 pouches per case
PSWE0002	PROSAT Sterile Polynit Heatseal LE Wipes Presaturated with 70% IPA and 30% WFI (EP), double bagged, flat stacked	300 x 300mm	30 wipes per pouch; 15 pouches per case
PSWE0003	PROSAT Sterile Polynit Heatseal LE Wipes Presaturated with 70% IPA and 30% WFI (EP), triple bagged	300 x 300mm	10 wipes per pouch; 30 pouches per case

For more information or to request a sample, email [infoeu@contecinc.com](mailto:infoeu@contecinc.com) or phone + 33 (0) 2 97 43 76 98.

Copyright© 2019 Contec, Inc. All rights reserved. INT042 071119

Contec  
Vannes Cedex  
France  
Tel: +33(0) 2 97 43 76 98

Contec Inc  
Spartanburg, SC  
United States  
Tel: +1 (864) 503-8333

Contec Cleanroom Technology (Suzhou) Co. Ltd  
Suzhou  
China  
Tel: +86-512-6274 4050

