

Contec® Sterile 70% Isopropyl Alcohol

USP <797> and USP <800> compliant sterile alcohol

Contec Sterile 70% Isopropanol contains 70% by volume USP Grade Isopropanol (IPA) and 30% USP Grade Purified Water. Filtered to a 0.2µ and induction sealed to assure sterility, the bottles are double bagged and gamma irradiated. Each bottle has a lot code and expiration date for easy record keeping. Product expires 5 years after date of manufacture. Validated per the Association for the Advancement of Medical Instrumentation (AAMI) 11137:2006 Guidelines to a 10⁻⁶ Sterility Assurance Level (SAL). Note: For spray configurations, Sterility Validation is performed after trigger is installed.

A convenient trigger sprayer with settings for stream or mist spraying is included with 16 and 32 oz. bottles where noted. The stream setting is ideal for targeting specific hard to reach areas, while the mist setting is used to cover larger areas. Each bottle is induction sealed to eliminate the occurrence of leakage. Instructions for sterile installation of the trigger sprayer are available. Both 16 oz. and 32 oz. USP <800> compliant flip top configurations are now available.



Use for wipe down of isolators, cabinets, laminar flow hoods and surfaces.

Features	Benefits
Contains 70% USP grade isopropyl alcohol	<ul style="list-style-type: none"> • Compliant with USP <797> and USP <800> • Evaporates and leaves extremely low residues • Eliminates mixing, processing and QC expense of in-house production of sterile isopropyl alcohol water blends
Filtered through a 0.2 micron in-line filter	<ul style="list-style-type: none"> • Eliminates the need for in-house filtration • Meets the requirement of USP grade purified water
Gamma irradiated to a Sterility Assurance Level (SAL) of 10 ⁻⁶	<ul style="list-style-type: none"> • Compliant with USP <797> and USP <800> • Packaged and sterilized for easy introduction into critical environments • Certificate of processing and sterility with each case
Induction sealed bottles	<ul style="list-style-type: none"> • Reduces opportunity for leakage during sterilization and shipment
Adjustable trigger sprayer	<ul style="list-style-type: none"> • Select spray or jet option
Flip-top cap	<ul style="list-style-type: none"> • Flip top for applications where spraying is not recommended
Plastic bottle	<ul style="list-style-type: none"> • Recyclable, environmentally friendly, no costly disposal
Double bagged	<ul style="list-style-type: none"> • To preserve product integrity
Lot numbered bottles	<ul style="list-style-type: none"> • Ensures complete product traceability

Part no.	Description		Size	Packaging
SB167030IR	Contec Sterile 70% Isopropyl alcohol with trigger sprayer	 Sterile	16 oz. (0.47 L)	12/case
HCFT7030IR	Contec Sterile 70% Isopropyl alcohol with flip-top cap	 Sterile	16 oz. (0.47 L)	12/case
HCFT7030IR-32	Contec Sterile 70% Isopropyl alcohol with flip-top cap	 Sterile	32 oz. (0.9 L)	12/case
SB327030IR	Contec Sterile 70% Isopropyl alcohol with trigger sprayer	 Sterile	32 oz. (0.9 L)	12/case
SB1287030IR	Contec Sterile 70% Isopropyl alcohol	 Sterile	1 gallon (3.8 L)	4/case