



# Irradiated Bucket Liners

## Description

BKTBAG2 bucket liners fit Contec's 15-liter polypropylene buckets and 10-liter stainless steel buckets. Gamma irradiated at 25-50 kGy, they are suitable for use in critical environments, including aseptic areas. Intended for single use, irradiated bucket liners minimize the risk of cross contamination and degradation of disinfectant solutions. The liners are constructed of cleanroom grade low-density polyethylene (LLDPE) film and sized 26" x 18" (66x46cm).



Irradiated Bucket Liners, fits 3.9-gal (15L) polypropylene buckets and 2.6-gal (10L) stainless steel buckets, double-bagged

Part Number	Description	EA/OB1	OB1/OB2	OB2/CS	EA/CS
BKTBAG2	Irradiated Bucket Liners, fits 3.9-gal (15L) polypropylene Buckets, double bagged	3	1	20	60

### Notes:

- 1) The information presented here is applicable to the part numbers shown below, as well as to any product containing the same materials and produced under the same conditions, regardless of product size or packaging configuration. Please contact a Contec sales representative for more details.
- 2) Data shown are typical values and should not be used as product specifications.
- 3) Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- 4) Current and/or comparison data may be available. Please contact a Contec sales representative for details.

Revision date: 03/02/18



Copyright © 2018 Contec, Inc. All rights reserved.

**North America**  
P.O. Box 530  
Spartanburg, SC 29304 USA  
Phone: +1-864-503-8333  
www.contecinc.com

**Europe**  
ZI du Prat – RP 3707  
56037 VANNES Cedex  
FRANCE  
+33 297 437 698

**China**  
No.17 Longyun Road  
Suzhou 215024  
CHINA  
+ 86-512-6274 4050